

## Fits your wheelchair the kit


The kit fits on different kinds of wheelchairs whether you have round or oval framework on your chair. 4 brackets mounted on the wheelchair and they are clamped around the frame. There will be no damage to the wheelchair.

### Wheelchair information:

I need to know how wide the front framework is in your wheelchair, so measure the length of 1A and 2A (see the picture). Measure from the center of each tube in millimeters.

The diameter of your frame in millimetres where the mounts 1 and 2 will be attached. If the tubes are oval shape, measure both the narrow and the wide side.

The distance between the brackets 1 and 2 can be between: 50-80mm.



### There are 3 sizes of ambike kit

Small: 1A and 2A are between 24-34 cm.  
Medium: 1A and 2A are between 34-45 cm.  
Large: 1A and 2A are between 45-54 cm.  
The maximum weight drivers 95kg

### Wheelchair information

The wheelchair's width 1A.....  
The wheelchair's width 2A.....  
The diameter of your frame 1.....  
The diameter of your frame 2.....  
If the tubes are oval shape, measure both the narrow and the wide side.

### Shipping Info Needed:

Full delivery addresses and phone number.  
The shipping Company needs this information for accurate delivery.  
You will receive an order confirmation form when all your info is received.  
Send all information to: [arm.bike@hotmail.com](mailto:arm.bike@hotmail.com)

**Do you change wheelchair and keeps you in the same measure as previously type Medium which is adjustable from 34 to 45cm. You do not change the frame. Change diameter of the wheelchair frame, you have to replace the brackets 1 and 2.**